

APPROVED	O.G. FIG.
BY	CLASS SUBCLASS
DRAFTSMAN	

Figure 1a. A general mixture synthesis with fluororous tags using a mixture of tagged compounds.




Figure 1b. A general mixture synthesis with fluororous tags using a mixture of tagged compounds and a mixture of reactants.






Figure 1c. A general mixture synthesis with fluororous tags using fluororous tagged reactants and a substrate.



APPROVED	O. G. FIG.	
BY	CLASS	SUBCLASS
DRAFTSMAN		


Figure 2. A representative example of a synthesis with a mixture of fluorine tagged compounds and a mixture of reactants



					Library	Esters	Products in order of retention times (min) on Fluorix column			
Tags		Esters		Thiols	1	1x, 2y, 3z	1xc (18.5); 1xb (18.9); 1xa (19.3); 1xd (23.8); 2yc (28.7); 2yb (28.7); 2ya (29.5); 2yd (32.6); 3zc (34.1); 3zb (34.1); 3za (35.1); 3zd (37.9)			
R ^f		R ¹	R ²							
x	C ₆ F ₁₃	1	Me	H	2	1y, 2z, 3x	3xc (18.1); 3xb (18.5); 3xa (18.7); 3xd (23.4); 1yc (27.0); 1yb (27.0); 1yc (27.6); 1yd (31.2); 2zc (35.6); 2zb (35.6); 2za (36.5); 2zd (38.8)			
y	C ₈ F ₁₇	2	Pr	H						
y	C ₁₀ F ₂₁	3	H	Me						
z				a C ₆ H ₅ b 2-naphthyl c p-MeOC ₆ H ₄ d p- ^t BuC ₆ H ₄						
					3	1z, 2x, 3y	2xc (20.4); 2xb (20.9); 2xa (21.0); 2xd (25.3); 3yc (26.4); 3yb (26.4); 3ya (27.0); 3yd (30.8); 1zc (34.2); 1zb (34.2); 1za (35.1); 1zd (37.8)			

APPROVED	O. G. FIG.
BY	CLASS SUBCLASS
DRAFTSMAN	

Figure 3. A representative HPLC trace of a library of compounds produced in the synthesis of Figure 2.^a



a) Retention times are listed in minutes; compound numbers refer to Figure 2; Fluofix column eluting with a gradient of 80% methanol/water increased to 100% methanol over 40 min. The peak at 39 min is an unknown impurity.

Figure 4. Preparation of Precursors for a Mixture Synthesis of Mappicine Analogs





Figure 5. Mappicine Mixture Synthesis and Separation



Mixture of **5a-e**

Mixture of protected mappicines **6a-e**

	R^1	R^2	R^3	R^4
a	SiMe_2^tBu	CH_2Ph	Me	C_6H_{13}
b	H	Et	ipr	C_4F_9
c	H	^tBu	ipr	C_6F_{13}
d	H	CH_2Ph	ipr	C_8F_{17}
e	Ph	Et	ipr	$\text{C}_{10}\text{F}_{21}$

separate on Fluofix™

Individual, pure samples of **6a,e**

Time	Gradient
0-5 min	80% MeOH/ H_2O -
5-25 min	90% MeOH/ H_2O
>25 min	100% MeOH

	Retention Time	Yield
a	3 min	36%
b	13 min	41%
c	18 min	29%
d	21 min	36%
e	28 min	43%

Figure 6. Preparation of precursors for a mixture synthesis of mappicine analogs (Example 8)

